

OCULUS BIOM[®] Optic Set

High-definition wide-field front lens and reduction lens in single-use design




The perfect view in every case

- ✓ Extremely wide field of view with outstanding resolution
- ✓ HD clarity under high magnification reduces the need for a contact lens
- ✓ Full-field clarity for decreased scleral depressing and panretinal laser
- ✓ Excellent depth of field for better stereopsis
- ✓ Improved view during air/fluid exchanges

The NEW BIOM[®] Optic Set

- Always ready, convenient blister pack design
- Minimizes risk of infection and cross contamination
- Same superior HD quality for each operation




Convenient and sterile blister package*



BIOM[®] HD Disposable Lens



New reduction lens increases surgical efficiency



Technical Data

Objective lens	BIOM [®]	BIOM [®] Optic Set (disposable)
f = 175 mm	BIOM [®] 4 m/c, 5 m/c	HD Professional for f = 175 mm
f = 200 mm	BIOM [®] 4 ml/cl, 5 ml/cl	HD Professional for f = 200 mm
	BIOM [®] 4 m/c, 5 m/c	HD Flex

* Blister contents: single-use reduction lens, thumbscrew cover, two drive belts for motorized BIOM[®], BIOM[®] HD Disposable Lens.



WWW.OCULUS.DE



The OCULUS QM system is certified in accordance with ISO 13485 (MDSAP) and (EU) 2017/745 (MDR)

OCULUS Optikgeräte GmbH
Postfach • 35549 Wetzlar • GERMANY
Tel. +49 641 2005-0 • Fax +49 641 2005-295
Email: export@oculus.de • www.oculus.de

Find your local OCULUS representative on our website.